


Product datasheet for **GA115373**

Human S100A16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100A16 CRISPRa kit - CRISPR gene activation of human S100 calcium binding protein A16
Symbol:	S100A16
UniProt ID:	Q96FQ6
Synonyms:	AAG13; DT1P1A7; S100F
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)