

Product datasheet for **GA115377**

Human Selenoprotein M (SELENOM) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SELENOM CRISPRa kit - CRISPR gene activation of human selenoprotein M
Symbol:	SELENOM
Components:	GA115377G1, Selenoprotein M gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115377G2, Selenoprotein M gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115377G3, Selenoprotein M gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WWX9
Synonyms:	SEPM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)