

Product datasheet for **GA115383**

Human C20orf96 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C20orf96 CRISPRa kit – CRISPR gene activation of human chromosome 20 open reading frame 96
Symbol:	C20orf96
Components:	GA115383G1, C20orf96 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115383G2, C20orf96 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115383G3, C20orf96 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NUD7</u>
Synonyms:	dJ1103G7.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

