

Product datasheet for **GA115386**

Human WFDC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WFDC3 CRISPRa kit - CRISPR gene activation of human WAP four-disulfide core domain 3
Symbol:	WFDC3
UniProt ID:	Q8IUW2
Synonyms:	dj447f3.3; WAP14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)