

Product datasheet for **GA115393**

Human **SAMD10** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SAMD10 CRISPRa kit - CRISPR gene activation of human sterile alpha motif domain containing 10
Symbol:	SAMD10
Components:	GA115393G1, SAMD10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115393G2, SAMD10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115393G3, SAMD10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BYL1</u>
Synonyms:	C20orf136; dj591C20; dj591C20.7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)