

Product datasheet for **GA115397**

Human C20orf117 (SOGA1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SOGA1 CRISPRa kit - CRISPR gene activation of human suppressor of glucose, autophagy associated 1
Symbol:	SOGA1
Components:	GA115397G1, C20orf117 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115397G2, C20orf117 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115397G3, C20orf117 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q94964</u>
Synonyms:	C20orf117; KIAA0889; SOGA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)