


Product datasheet for **GA115398**

Human C20orf118 (TLDC2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TLDC2 CRISPRa kit - CRISPR gene activation of human TBC/LysM-associated domain containing 2
Symbol:	TLDC2
Components:	GA115398G1, C20orf118 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115398G2, C20orf118 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115398G3, C20orf118 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A0PIX2
Synonyms:	C20orf118
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)