

## Product datasheet for **GA115404**

### Human TMEM37 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | TMEM37 CRISPRa kit - CRISPR gene activation of human transmembrane protein 37                                                                                                                                                                                                                                 |
| Symbol:       | TMEM37                                                                                                                                                                                                                                                                                                        |
| Components:   | GA115404G1, TMEM37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115404G2, TMEM37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115404G3, TMEM37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8WXS4</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | PR; PR1                                                                                                                                                                                                                                                                                                       |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)