


Product datasheet for **GA115424**

Human RIPOR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RIPOR3 CRISPRa kit - CRISPR gene activation of human RIPOR family member 3
Symbol:	RIPOR3
Components:	GA115424G1, RIPOR3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115424G2, RIPOR3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115424G3, RIPOR3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96MK2</u>
Synonyms:	C20orf175; C20orf176
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)