

Product datasheet for **GA115433**

Human **STK35** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	STK35 CRISPRa kit - CRISPR gene activation of human serine/threonine kinase 35
Symbol:	STK35
UniProt ID:	Q8TDR2
Synonyms:	CLIK1; STK35L1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)