

Product datasheet for **GA115439**

Human ITLN2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITLN2 CRISPRa kit - CRISPR gene activation of human intelectin 2
Symbol:	ITLN2
Components:	GA115439G1, ITLN2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115439G2, ITLN2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115439G3, ITLN2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WWU7</u>
Synonyms:	HL-2; HL2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)