


Product datasheet for **GA115440**

Human **RAB40A** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	RAB40A CRISPRa kit - CRISPR gene activation of human RAB40A, member RAS oncogene family
Symbol:	RAB40A
Components:	GA115440G1, RAB40A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115440G2, RAB40A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115440G3, RAB40A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WXH6</u>
Synonyms:	RAR2; RAR2A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)