


Product datasheet for **GA115458**

Human LDLRAD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LDLRAD3 CRISPRa kit - CRISPR gene activation of human low density lipoprotein receptor class A domain containing 3
Symbol:	LDLRAD3
Components:	GA115458G1, LDLRAD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115458G2, LDLRAD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115458G3, LDLRAD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86YD5</u>
Synonyms:	LRAD3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)