

Product datasheet for **GA115459**

Human PSMA8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMA8 CRISPRa kit - CRISPR gene activation of human proteasome subunit alpha 8
Symbol:	PSMA8
UniProt ID:	<u>Q8TAA3</u>
Synonyms:	PSMA7L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)