

## Product datasheet for **GA115468**

### Human **FAM76B** activation kit by **CRISPRa**

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | FAM76B CRISPRa kit - CRISPR gene activation of human family with sequence similarity 76 member B                                                                                                                                                                                                              |
| <b>Symbol:</b>       | FAM76B                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA115468G1, FAM76B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115468G2, FAM76B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115468G3, FAM76B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q5HYJ3</a></u>                                                                                                                                                                                                                                                                                 |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

