

Product datasheet for **GA115470**

Human PIWIL4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIWIL4 CRISPRa kit - CRISPR gene activation of human piwi like RNA-mediated gene silencing 4
Symbol:	PIWIL4
UniProt ID:	Q7Z3Z4
Synonyms:	HIWI2; MIWI2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)