

Product datasheet for **GA115474**

Human **POGLUT3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	POGLUT3 CRISPRa kit - CRISPR gene activation of human protein O-glucosyltransferase 3
Symbol:	POGLUT3
Components:	GA115474G1, POGLUT3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115474G2, POGLUT3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115474G3, POGLUT3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7Z4H8
Synonyms:	MGC33424
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

