

Product datasheet for **GA115477**

Human SPINDOC activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPINDOC CRISPRa kit - CRISPR gene activation of human spindlin interactor and repressor of chromatin binding
Symbol:	SPINDOC
Components:	GA115477G1, SPINDOC gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115477G2, SPINDOC gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115477G3, SPINDOC gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BUA3</u>
Synonyms:	DKFZp762N0114
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)