

Product datasheet for **GA115478**

Human PLEKHA7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA7 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A7
Symbol:	PLEKHA7
UniProt ID:	Q6IQ23
Synonyms:	DKFZp686M22243
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)