


Product datasheet for **GA115479**

Human SPTY2D1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPTY2D1 CRISPRa kit - CRISPR gene activation of human SPT2 chromatin protein domain containing 1
Symbol:	SPTY2D1
Components:	GA115479G1, SPTY2D1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115479G2, SPTY2D1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115479G3, SPTY2D1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68D10</u>
Synonyms:	Spt2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)