

Product datasheet for **GA115483**

Human DNHD1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DNHD1 CRISPRa kit - CRISPR gene activation of human dynein heavy chain domain 1
Symbol:	DNHD1
Components:	GA115483G1, DNHD1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115483G2, DNHD1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115483G3, DNHD1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96M86
Synonyms:	C11orf47; CCDC35; DHCD1; DNHD1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)