


Product datasheet for **GA115495**

Human CFAP251 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR66 CRISPRa kit - CRISPR gene activation of human WD repeat domain 66
Symbol:	CFAP251
Components:	GA115495G1, CFAP251 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115495G2, CFAP251 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115495G3, CFAP251 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TBY9</u>
Synonyms:	CaM-IP4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)