

Product datasheet for **GA115497**

Human TSPAN19 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPAN19 CRISPRa kit – CRISPR gene activation of human tetraspanin 19
Symbol:	TSPAN19
Components:	GA115497G1, TSPAN19 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115497G2, TSPAN19 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115497G3, TSPAN19 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P0C672
Synonyms:	FLJ44351
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

