

## Product datasheet for **GA115498**

### Human Bestrophin 3 (BEST3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	BEST3 CRISPRa kit - CRISPR gene activation of human bestrophin 3
Symbol:	BEST3
UniProt ID:	<u><a href="#">Q8N1M1</a></u>
Synonyms:	VMD2L3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)