

Product datasheet for **GA115515**

Human Proline rich membrane anchor 1 (PRIMA1) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PRIMA1 CRISPRa kit - CRISPR gene activation of human proline rich membrane anchor 1 |
| Symbol: | PRIMA1 |
| UniProt ID: | Q86XR5 |
| Synonyms: | acetylcholinesterase membrane anchor precursor PRiMA; OTTHUMP00000180666; PRIMA; proline rich membrane anchor 1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)