

Product datasheet for **GA115517**

Human C14orf50 (PPP1R36) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPP1R36 CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 36
Symbol:	PPP1R36
Components:	GA115517G1, C14orf50 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115517G2, C14orf50 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115517G3, C14orf50 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96LQ0</u>
Synonyms:	C14orf50
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)