

Product datasheet for **GA115523**

Human C14orf166B (LRRC74A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC74A CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 74A
Symbol:	LRRC74A
Components:	GA115523G1, C14orf166B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115523G2, C14orf166B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115523G3, C14orf166B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q0VAA2
Synonyms:	C14orf166B; LRRC74
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)