

## Product datasheet for **GA115537**

### Human HAPLN3 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | HAPLN3 CRISPRa kit - CRISPR gene activation of human hyaluronan and proteoglycan link protein 3                                                                                                                                                                                                               |
| Symbol:       | HAPLN3                                                                                                                                                                                                                                                                                                        |
| Components:   | GA115537G1, HAPLN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115537G2, HAPLN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA115537G3, HAPLN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q96S86</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | EXLD1; HsT19883                                                                                                                                                                                                                                                                                               |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)