

Product datasheet for **GA115539**

Human TMC05A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMC05A CRISPRa kit - CRISPR gene activation of human transmembrane and coiled-coil domains 5A
Symbol:	TMC05A
UniProt ID:	<u>Q8N6Q1</u>
Synonyms:	TMC05
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)