


Product datasheet for **GA115553**

Human TEK5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEK5 CRISPRa kit - CRISPR gene activation of human tektin 5
Symbol:	TEK5
Components:	GA115553G1, TEK5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115553G2, TEK5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115553G3, TEK5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96M29</u>
Synonyms:	CT149
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)