

Product datasheet for **GA115560**

Human Interleukin 34 (IL34) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL34 CRISPRa kit - CRISPR gene activation of human interleukin 34
Symbol:	IL34
UniProt ID:	<u>Q6ZMJ4</u>
Synonyms:	C16orf77; IL-34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)