

Product datasheet for **GA115563**

Human TMED6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMED6 CRISPRa kit - CRISPR gene activation of human transmembrane p24 trafficking protein 6
Symbol:	TMED6
Components:	GA115563G1, TMED6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115563G2, TMED6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115563G3, TMED6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WW62
Synonyms:	p24g5; PRO34237; SPL9146
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)