


Product datasheet for **GA115571**

Human TEPSIN activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TEPSIN CRISPRa kit - CRISPR gene activation of human TEPSIN adaptor related protein complex 4 accessory protein
Symbol:	TEPSIN
Components:	GA115571G1, TEPSIN gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115571G2, TEPSIN gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115571G3, TEPSIN gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96N21</u>
Synonyms:	C17orf56; ENTHD2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)