


Product datasheet for **GA115574**

Human CD300 (CD300LF) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | CD300LF CRISPRa kit - CRISPR gene activation of human CD300 molecule like family member f |
| Symbol: | CD300LF |
| UniProt ID: | <u>Q8TDQ1</u> |
| Synonyms: | CD300f; CLM-1; CLM1; IgSF13; IREM-1; IREM1; LMIR3; NKIR |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)