

Product datasheet for **GA115575**

Human CD300LD-AS1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C17orf77 CRISPRa kit - CRISPR gene activation of human chromosome 17 open reading frame 77
Symbol:	CD300LD-AS1
Components:	GA115575G1, CD300LD-AS1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115575G2, CD300LD-AS1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115575G3, CD300LD-AS1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96MU5</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)