


Product datasheet for **GA115583**

Human PIK3R6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3R6 CRISPRa kit - CRISPR gene activation of human phosphoinositide-3-kinase regulatory subunit 6
Symbol:	PIK3R6
Components:	GA115583G1, PIK3R6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115583G2, PIK3R6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115583G3, PIK3R6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5UE93</u>
Synonyms:	C17orf38; HsT41028; p84 PIKAP; p87(PIKAP); p87PIKAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)