

Product datasheet for **GA115593**

Human TMEM199 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM199 CRISPRa kit - CRISPR gene activation of human transmembrane protein 199
Symbol:	TMEM199
Components:	GA115593G1, TMEM199 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115593G2, TMEM199 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115593G3, TMEM199 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N511</u>
Synonyms:	C17orf32; CDG2P; VMA12; VPH2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)