

Product datasheet for **GA115596**

Human KCTD11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCTD11 CRISPRa kit - CRISPR gene activation of human potassium channel tetramerization domain containing 11
Symbol:	KCTD11
UniProt ID:	Q693B1
Synonyms:	C17orf36; KCASH1; REN; REN/KCTD11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)