

Product datasheet for **GA115600**

Human WIPF2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WIPF2 CRISPRa kit - CRISPR gene activation of human WAS/WASL interacting protein family member 2
Symbol:	WIPF2
Components:	GA115600G1, WIPF2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115600G2, WIPF2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115600G3, WIPF2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TF74</u>
Synonyms:	WICH; WIRE
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)