

Product datasheet for **GA115624**

Human KLC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KLC3 CRISPRa kit - CRISPR gene activation of human kinesin light chain 3
Symbol:	KLC3
Components:	GA115624G1, KLC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115624G2, KLC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115624G3, KLC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P597</u>
Synonyms:	kinesin light chain 2-like; kinesin light chain 3; KLC2; KLC2, KLCt, KLC2L, KNS2B; KLC2-like; KLC2L; KLCt; KNS2B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)