

Product datasheet for **GA115628**

Human TMEM190 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM190 CRISPRa kit - CRISPR gene activation of human transmembrane protein 190
Symbol:	TMEM190
Components:	GA115628G1, TMEM190 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115628G2, TMEM190 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115628G3, TMEM190 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WZ59</u>
Synonyms:	MDAC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)