

## Product datasheet for **GA115630**

### Human TMC4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMC4 CRISPRa kit - CRISPR gene activation of human transmembrane channel like 4
Symbol:	TMC4
UniProt ID:	<u><a href="#">Q7Z404</a></u>
Synonyms:	MGC39329
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)