

Product datasheet for **GA115639**

Human ZNF420 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF420 CRISPRa kit - CRISPR gene activation of human zinc finger protein 420
Symbol:	ZNF420
UniProt ID:	<u>Q8TAQ5</u>
Synonyms:	APAK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)