

Product datasheet for **GA115646**

Human DPY19L3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DPY19L3 CRISPRa kit - CRISPR gene activation of human dpy-19 like C-mannosyltransferase 3
Symbol:	DPY19L3
Components:	GA115646G1, DPY19L3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115646G2, DPY19L3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115646G3, DPY19L3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZPD9</u>
Synonyms:	DKFZp686j17135; MGC35440
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)