

Product datasheet for **GA115662**

Human ATP8B3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ATP8B3 CRISPRa kit - CRISPR gene activation of human ATPase phospholipid transporting 8B3
Symbol:	ATP8B3
Components:	GA115662G1, ATP8B3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115662G2, ATP8B3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115662G3, ATP8B3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O60423</u>
Synonyms:	ATPIK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)