

Product datasheet for **GA115675**

Human PHF13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHF13 CRISPRa kit - CRISPR gene activation of human PHD finger protein 13
Symbol:	PHF13
Components:	GA115675G1, PHF13 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115675G2, PHF13 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115675G3, PHF13 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86YI8</u>
Synonyms:	PHF5; SPOC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)