

Product datasheet for **GA115691**

Human MOBKL2C (MOB3C) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MOB3C CRISPRa kit - CRISPR gene activation of human MOB kinase activator 3C
Symbol:	MOB3C
UniProt ID:	<u>Q70IA8</u>
Synonyms:	MOB1E; MOBKL2C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)