


Product datasheet for **GA115694**

Human LELP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LELP1 CRISPRa kit - CRISPR gene activation of human late cornified envelope like proline rich 1
Symbol:	LELP1
Components:	GA115694G1, LELP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115694G2, LELP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115694G3, LELP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T871</u>
Synonyms:	late cornified envelope-like proline-rich 1; novel small proline rich protein; OTTHUMP00000015821
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)