

## Product datasheet for **GA115701**

### Human IL23 Receptor (IL23R) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	IL23R CRISPRa kit - CRISPR gene activation of human interleukin 23 receptor
Symbol:	IL23R
UniProt ID:	<u><a href="#">Q5VWK5</a></u>
Synonyms:	interleukin 23 receptor; OTTHUMP00000010660
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)