

## Product datasheet for **GA115704**

### Human SHISA4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SHISA4 CRISPRa kit - CRISPR gene activation of human shisa family member 4
Symbol:	SHISA4
Components:	GA115704G1, SHISA4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA115704G2, SHISA4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA115704G3, SHISA4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q96DD7</a></u>
Synonyms:	C1orf40; TMEM58
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)